

Atypical Chemokine Receptor D6 protein (STJ503623)

STJ503623

GENERAL INFORMATION

Product Type	Proteins
Short Description	Atypical Chemokine Receptor D6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59289
Gene Symbol	Ackr2
Uniprot ID	ACKR2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 300-378 mouse Atypical chemokine receptor 2 protein.
Immunogen Region	300-378
Specificity	
Immunogen Sequence	