

Alpha 1 Fetoprotein protein (STJ503631)

STJ503631

GENERAL INFORMATION

Product Type	Proteins
Short Description	Alpha 1 Fetoprotein-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	174
Gene Symbol	AFP
Uniprot ID	FETA_HUMAN
Immunogen	Full length native protein (purified) corresponding to Human alpha 1 Fetoprotein. Purified human alpha 1 Fetoprotein.
Immunogen Region	Full length
Specificity	
Immunogen Sequence	