

AIM2 protein (STJ503633)
STJ503633**GENERAL INFORMATION**

Product Type	Proteins
Short Description	AIM2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9447
Gene Symbol	AIM2
Uniprot ID	AIM2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 240-340 on human AIM2 protein.
Immunogen Region	240-340
Specificity	
Immunogen Sequence	