

AKR1B1 protein (STJ503635)

STJ503635

GENERAL INFORMATION

Product Type	Proteins
Short Description	AKR1B1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	231
Gene Symbol	AKR1B1
Uniprot ID	ALDR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-350 on human Aldo-keto reductase family 1, member B1 protein.
Immunogen Region	250-350
Specificity	
Immunogen Sequence	