

AMACR protein (STJ503641)

STJ503641

GENERAL INFORMATION

Product Type	Proteins
Short Description	AMACR-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100 ICC: 1:100 IF: 1:200 CFM: 1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23600
Gene Symbol	AMACR
Uniprot ID	AMACR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human alpha-methylacyl-CoA racemase protein.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	