

Annexin A10 protein (STJ503650)

STJ503650

GENERAL INFORMATION

Product Type	Proteins
Short Description	Annexin A10-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11199
Gene Symbol	ANXA10
Uniprot ID	ANX10_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-100 on human Annexin A10 protein.
Immunogen Region	60-100
Specificity	
Immunogen Sequence	