

Annexin A6 protein (STJ503652)

STJ503652

GENERAL INFORMATION

Product Type	Proteins
Short Description	Annexin A6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	309
Gene Symbol	ANXA6
Uniprot ID	ANXA6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 500-560 on human Annexin A6 protein.
Immunogen Region	500-560
Specificity	
Immunogen Sequence	