

Aph 1b protein (STJ503656)

STJ503656

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aph 1b-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83464
Gene Symbol	APH1B
Uniprot ID	APH1B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 235-257 of Human Presenilin stabilizing protein.
Immunogen Region	235-257
Specificity	
Immunogen Sequence	