

Aquaporin 11 protein (STJ503660)

STJ503660

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aquaporin 11-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	286758
Gene Symbol	Aqp11
Uniprot ID	AQP11_RAT
Immunogen	Synthetic peptide taken within amino acid region 140-190 on Aquaporin 11 protein.
Immunogen Region	140-190
Specificity	
Immunogen Sequence	