

Aquaporin 6 protein (STJ503666)

STJ503666

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aquaporin 6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29170
Gene Symbol	Aqp6
Uniprot ID	AQP6_RAT
Immunogen	Synthetic peptide taken within amino acid region 220-270 on rat Aquaporin 6 protein.
Immunogen Region	220-270
Specificity	
Immunogen Sequence	