

Aquaporin 7 protein (STJ503667)

STJ503667

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aquaporin 7-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11832
Gene Symbol	Aqp7
Uniprot ID	AQP7_MOUSE
Immunogen	Synthetic amidated peptide taken within amino acid region 150-200 on rat Aquaporin 7 protein, covalently modified and conjugated to KLH.
Immunogen Region	150-200
Specificity	
Immunogen	
Sequence	