

Aquaporin 8 protein (STJ503668)

STJ503668

GENERAL INFORMATION

Product Type	Proteins
Short Description	Aquaporin 8-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29172
Gene Symbol	Aqp8
Uniprot ID	AQP8_RAT
Immunogen	Synthetic peptide taken within amino acid region 215-263 on rat Aquaporin 8 protein.
Immunogen Region	215-263
Specificity	
Immunogen Sequence	