

ASF1A protein (STJ503686)
STJ503686**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ASF1A-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25842
Gene Symbol	ASF1A
Uniprot ID	ASF1A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-120 on human AF1A protein.
Immunogen Region	60-120
Specificity	
Immunogen Sequence	