

Fbx32 protein (STJ503688)
STJ503688**GENERAL INFORMATION**

Product Type	Proteins
Short Description	Fbx32-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114907
Gene Symbol	FBXO32
Uniprot ID	FBX32_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 240-280 on Atr1 protein.
Immunogen Region	240-280
Specificity	
Immunogen Sequence	