

Axin 2 protein (STJ503690)
STJ503690**GENERAL INFORMATION**

Product Type	Proteins
Short Description	Axin 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8313
Gene Symbol	AXIN2
Uniprot ID	AXIN2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 800-843 on human Axin 2 protein.
Immunogen Region	800-843
Specificity	
Immunogen Sequence	