

Axl protein (STJ503691)

STJ503691

GENERAL INFORMATION

Product Type	Proteins
Short Description	Axl-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26362
Gene Symbol	Axl
Uniprot ID	UFO_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 840-880 on mouse AXL protein. The peptide was post translationally modified to achieve proper orientation and antigenicity.
Immunogen Region	840-880
Specificity	
Immunogen Sequence	