

BCKDH protein (STJ503694)

STJ503694

GENERAL INFORMATION

Product Type	Proteins
Short Description	BCKDH-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10295
Gene Symbol	BCKDK
Uniprot ID	BCKD_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 170-200 on human BCKD protein.
Immunogen Region	170-200
Specificity	
Immunogen Sequence	