

Bcl6 protein (STJ503695)
STJ503695**GENERAL INFORMATION**

Product Type	Proteins
Short Description	Bcl6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	604
Gene Symbol	BCL6
Uniprot ID	BCL6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 610-660 on human B-cell lymphoma protein isoform 1.
Immunogen Region	610-660
Specificity	
Immunogen Sequence	