

SLC7A10 protein (STJ503698)

STJ503698

GENERAL INFORMATION

Product Type	Proteins
Short Description	SLC7A10-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56301
Gene Symbol	SLC7A10
Uniprot ID	AAA1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-500 on rat Asc-type amino acid transporter 1 protein.
Immunogen Region	450-500
Specificity	
Immunogen Sequence	