

ASF1A protein (STJ503703)
STJ503703**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ASF1A-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25842
Gene Symbol	ASF1A
Uniprot ID	ASF1A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	