

Baboon protein (STJ503714)

STJ503714

GENERAL INFORMATION

Product Type	Proteins
Short Description	Baboon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	N-epitope Synthetic peptide corresponding to amino acid region 200-221 of Drosophila melanogaster protein. The peptide was covalently modified for conjugation. Synthetic peptide: CNKNENMKRIFEDYVPHRLTS
Immunogen	200-221
Region	
Specificity	
Immunogen	
Sequence	