

Baboon protein (STJ503716)

STJ503716

GENERAL INFORMATION

Product Type	Proteins
Short Description	Baboon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	C-epitope Synthetic peptide corresponding to amino acids 552-572 of Drosophila melanogaster Baboon protein. Peptide was covalently modified. Sequence: QPDPSIEEMKKVVCIEKCRPN
Immunogen Region	552-572
Specificity	
Immunogen Sequence	