

Bestrophin 2 protein (STJ503720)

STJ503720

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | Bestrophin 2-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|-------------------------------|---|
| Gene ID | 212989 |
| Gene Symbol | Best2 |
| Uniprot ID | BEST2_MOUSE |
| Immunogen | Synthetic peptide taken within amino acid region 388-404 of mouse Bestrophin 2 protein, peptide was covalently modified. Sequence: EVHGDFLQRLLPAGAGS |
| Immunogen Region | 388-404 |
| Specificity | |
| Immunogen Sequence | |