

Bestrophin 3 protein (STJ503722)

STJ503722

GENERAL INFORMATION

Product Type	Proteins
Short Description	Bestrophin 3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	144453
Gene Symbol	BEST3
Uniprot ID	BEST3_HUMAN
Immunogen	Synthetic peptide taken within the last 50 amino acids from the C-terminus of human Bestrophin-3 protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	