

BEST4 protein (STJ503724)
STJ503724**GENERAL INFORMATION**

Product Type	Proteins
Short Description	BEST4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	266675
Gene Symbol	BEST4
Uniprot ID	BEST4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-500 on human Bestrophin 4 protein.
Immunogen Region	400-500
Specificity	
Immunogen Sequence	