

## BMP2 protein (STJ503727)

STJ503727

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	BMP2-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	12156
<b>Gene Symbol</b>	Bmp2
<b>Uniprot ID</b>	BMP2_MOUSE
<b>Immunogen</b>	Synthetic peptide corresponding to Mouse BMP2 amino acid region 45-60 conjugated to Keyhole Limpet Haemocyanin (KLH). Sequence: SEDVLSSEFELRLLSMF
<b>Immunogen</b>	45-60
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	