

BMP4 protein (STJ503728)

STJ503728

GENERAL INFORMATION

Product Type	Proteins
Short Description	BMP4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 184-197 on rat BMP4 protein, peptide was covalently modified. Sequence: VMKPPAPGHLIC
Immunogen Region	184-197
Specificity	
Immunogen Sequence	