

BMP8 protein (STJ503729)
STJ503729**GENERAL INFORMATION**

Product Type	Proteins
Short Description	BMP8-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	656
Gene Symbol	BMP8B
Uniprot ID	BMP8B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 160-200 on human BMP8 protein
Immunogen Region	160-200
Specificity	
Immunogen Sequence	