

BMP11 protein (STJ503730)
STJ503730**GENERAL INFORMATION**

Product Type	Proteins
Short Description	BMP11-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10220
Gene Symbol	GDF11
Uniprot ID	GDF11_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 60-100 on human BMP11/GDF11 protein.
Immunogen Region	60-100
Specificity	
Immunogen Sequence	