

Bone marrow stromal cell antigen 1 protein (STJ503731)

STJ503731

GENERAL INFORMATION

Product Type	Proteins
Short Description	Bone marrow stromal cell antigen 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	683
Gene Symbol	BST1
Uniprot ID	BST1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids 255-274 of Human Bone marrow stromal cell antigen 1. Peptide was covalently modified. Sequence: FQY SCI NDY RPV KLL QCV DH
Immunogen Region	255-274
Specificity	
Immunogen Sequence	