

BST2 protein (STJ503734)
STJ503734**GENERAL INFORMATION**

Product Type	Proteins
Short Description	BST2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	684
Gene Symbol	BST2
Uniprot ID	BST2_HUMAN
Immunogen	C-epitope Synthetic peptide taken within amino acid region 100-150 on human Bonemarrow stromal antigen 2 protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	