

BTK protein (STJ503737)

STJ503737

GENERAL INFORMATION

Product Type	Proteins
Short Description	BTK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	695
Gene Symbol	BTK
Uniprot ID	BTK_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 550-650 on human BTK protein
Immunogen Region	550-650
Specificity	
Immunogen Sequence	