

Beta crystallin S protein (STJ503739)

STJ503739

GENERAL INFORMATION

Product Type	Proteins
Short Description	Beta crystallin S-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12970
Gene Symbol	Crygs
Uniprot ID	CRYGS_MOUSE
Immunogen	C-epitope Synthetic peptide taken within amino acid region 120-178 on mouse BXTIS protein.
Immunogen Region	120-178
Specificity	
Immunogen Sequence	