

C1orf146 protein (STJ503756)

STJ503756

GENERAL INFORMATION

Product Type	Proteins
Short Description	C1orf146-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	388649
Gene Symbol	SPO16
Uniprot ID	SPO16_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-180 on human C1orf146 protein.
Immunogen Region	100-180
Specificity	
Immunogen Sequence	