

C2orf88 protein (STJ503758)

STJ503758

GENERAL INFORMATION

Product Type	Proteins
Short Description	C2orf88-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84281
Gene Symbol	C2orf88
Uniprot ID	SMAKA_HUMAN
Immunogen	Mid Region epitope synthetic peptide taken within amino acid region 25-100 on human C2orf88 protein.
Immunogen Region	25-100
Specificity	
Immunogen Sequence	