

C3Orf34 protein (STJ503760)

STJ503760

GENERAL INFORMATION

Product Type	Proteins
Short Description	C3Orf34-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84984
Gene Symbol	CEP19
Uniprot ID	CEP19_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 120-170 on human Chromosome open reading frame 34 protein.
Immunogen Region	120-170
Specificity	
Immunogen Sequence	