

C4a protein (STJ503761)

STJ503761

GENERAL INFORMATION

Product Type	Proteins
Short Description	C4a-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	720/721
Gene Symbol	C4A
Uniprot ID	CO4A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1220-1270 on human complement C4-A preprotein.
Immunogen Region	1220-1270
Specificity	
Immunogen Sequence	