

CaMKII gamma protein (STJ503771)

STJ503771

GENERAL INFORMATION

Product Type	Proteins
Short Description	CaMKII gamma-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	818
Gene Symbol	CAMK2G
Uniprot ID	KCC2G_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-455 on human CaMKII gamma protein.
Immunogen Region	400-455
Specificity	
Immunogen Sequence	