

CACNB2 protein (STJ503778)

STJ503778

GENERAL INFORMATION

Product Type	Proteins
Short Description	CACNB2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116600
Gene Symbol	Cacnb2
Uniprot ID	CACB2_RAT
Immunogen	Synthetic peptide taken within amino acid region 500-600 on rat voltage-dependent L-type calcium channel subunit beta-2 protein
Immunogen Region	500-600
Specificity	
Immunogen	
Sequence	