

CCK4 protein (STJ503779)

STJ503779

GENERAL INFORMATION

Product Type	Proteins
Short Description	CCK4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7786
Gene Symbol	MAP3K12
Uniprot ID	M3K12_HUMAN
Immunogen	Purified Protein tyrosine kinase 7 protein was emulsified in complete adjuvant and used to immunize the rabbits.
Immunogen Region	
Specificity	
Immunogen Sequence	