

CCR4 protein (STJ503784)

STJ503784

GENERAL INFORMATION

Product Type	Proteins
Short Description	CCR4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1233
Gene Symbol	CCR4
Uniprot ID	CCR4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C-C chemokine receptor type 4 protein.
Immunogen Region	1-50aa
Specificity	
Immunogen Sequence	