

CCR5 protein (STJ503785)

STJ503785

GENERAL INFORMATION

Product Type	Proteins
Short Description	CCR5-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1234
Gene Symbol	CCR5
Uniprot ID	CCR5_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human CCR5 protein.
Immunogen Region	1-50aa
Specificity	
Immunogen Sequence	