

CCR6 protein (STJ503786)
STJ503786**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CCR6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1235
Gene Symbol	CCR6
Uniprot ID	CCR6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human C-C chemokine receptor type 6 protein.
Immunogen Region	1-50aa
Specificity	
Immunogen Sequence	