

CD19 protein (STJ503791)

STJ503791

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD19-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	930
Gene Symbol	CD19
Uniprot ID	CD19_HUMAN
Immunogen	Synthetic peptide within Human CD19 amino acid region 228-242. The selected peptide was post-synthetically modified to achieve highest antigenicity before coupling to carrier protein using hetero bifunctional cross linker. Sequence: LELKDDRPARDMWMM
Immunogen	228-242
Region	
Specificity	
Immunogen	
Sequence	