

CD2 protein (STJ503793)

STJ503793

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	914
Gene Symbol	CD2
Uniprot ID	CD2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 3335-351 on human T-cell surface antigen CD2 protein.
Immunogen Region	3335-351
Specificity	
Immunogen Sequence	