

CD21 protein (STJ503794)
STJ503794**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CD21-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1380
Gene Symbol	CR2
Uniprot ID	CR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 473-523 on mouse GLUT2 protein.
Immunogen Region	473-523
Specificity	
Immunogen Sequence	