

## CD21 protein (STJ503795)

STJ503795

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	CD21-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	1380
<b>Gene Symbol</b>	CR2
<b>Uniprot ID</b>	CR2_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 915-965 on human Complement Receptor type 2 isoform X1 protein.
<b>Immunogen Region</b>	915-965
<b>Specificity</b>	
<b>Immunogen Sequence</b>	