

CD3 epsilon protein (STJ503798)

STJ503798

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD3 epsilon-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	916
Gene Symbol	CD3E
Uniprot ID	CD3E_HUMAN
Immunogen	C-epitope synthetic peptide corresponding to unique amino acid sequence on CD3 epsilon protein.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	