

CD3 Gamma protein (STJ503800)

STJ503800

GENERAL INFORMATION

Product Type	Proteins
Short Description	CD3 Gamma-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	917
Gene Symbol	CD3G
Uniprot ID	CD3G_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 131-181 on human T-cell surface glycoprotein CD3 gamma protein.
Immunogen Region	131-181
Specificity	
Immunogen Sequence	