

C1orf146 protein (STJ503801)

STJ503801

GENERAL INFORMATION

Product Type	Proteins
Short Description	C1orf146-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	388649
Gene Symbol	SPO16
Uniprot ID	SPO16_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	